G-IVF™ PLUS

26078.02

Innovative Cell and Tissue Technology

EN: Medium for preparation and handling of gametes and for in vitro fertilisation.

SUPPLEMENTED WITH HSA

G-IVF™ PLUS is a bicarbonate buffered medium containing human serum albumin and gentamicin as an antibacterial agent.

Ready to use after equilibration at +37 °C and 6 % CO₂ atmosphere.

Storage instructions and stability

Store dark at +2 to +8 °C.

G-IVF™ PLUS is stable until the expiry date shown on the container labels and the LOT-specific Certificate of Analysis.

Media bottles should not be stored after opening. Discard excess media after completion of the procedure.

Directions for use

The following is the general procedure for using G-IVF™ PLUS. For more information on the recommended use of G-IVF™ PLUS, please see "Vitrolife G Series™ Manual". The Manual is available on request from Vitrolife, and may be obtained through the company's internet website, www.vitrolife.com. Directions for use in other languages can also be found at the website.

Sperm preparation:

Swim-up migration method for semen samples with a good sperm count and motility.

Pipette 1.0 mL of semen into the tube. Overlay 2.0 mL of equilibrated G-IVF™ PLUS. Place the swim-up tube in an angled position in the incubator at +37 °C and 6 % CO₂ for 30–60 minutes. Aspirate the top medium without touching the underlying semen and transfer to a clean tube. Add 5.0 mL of G-IVF™ PLUS, mix and centrifuge for 10 minutes at 300–600g.

Aspirate the supernatant and re-suspend the pellet in 5.0 mL of equilibrated G-IVF™ PLUS and repeat the centrifugation. Discard the supernatant and resuspend the pellet in 0.5–1.0 mL of equilibrated G-IVF™ PLUS. Determine motility and concentration of spermatozoa in the washed sample.

Density gradient centrifugation method.

This method can be used to wash all samples regardless of quality.

Mix SpermGrad™ with G-IVF™ PLUS in separate tubes to obtain 90% and 45% stock solutions. For 90% stock solution, mix 9.0 mL SpermGrad™ with 1.0 mL G-IVF™ PLUS. For 45% stock solution mix 4.5 mL SpermGrad™ with 5.5 mL G-IVF™ PLUS. Mix the solutions thoroughly. Store in sterile non-toxic tubes or sterile tissue culture flasks. Label and refrigerate until use. Equilibrate the solutions at +37 °C and 6% CO₂ before use. The density gradient should be layered in a sterile and rinsed conical non-toxic centrifuge tube. Pipette 1.5

mL of 90 % solution into the tube first and then slowly 1.5 mL of 45 % solution on top of it. Finally, the semen is gently layered on the top. If the semen sample is of poor quality, up to 2 mL of semen can be layered on top of the gradient. Adding too much semen will result in poor separation. The tube is then centrifuged for 10–20 minutes at 300–600g. Remove the two top layers and take care not to leave any residues on the tube wall. Transfer the sperm pellet with as little of the 90 % solution as possible to a sterile conical rinsed tube with 5 mL of equilibrated G-IVF™ PLUS. Centrifuge for 10 minutes at 300–600g. Discard the supernatant and repeat the wash. After the second wash, the pellet is re-suspended in 1 mL of equilibrated G-IVF™ PLUS. The washed sample is then assessed for motility and concentration. Dilute the washed sample with equilibrated G-IVF™ PLUS to a final concentration of 75 000–200 000 motile sperm/mL. Prepare rinsed insemination dishes with 0.5–1.0 mL of sperm suspension and equilibrate at +37 °C and 6 % CO₂ for at least 2 hours.

Oocyte identification:

Rinse the oocytes first in pre-warmed supplemented G-MOPS™/G-MOPS™ PLUS, and then in equilibrated G-IVF™ PLUS. The washing procedure should include at least two steps with 1.0 mL of G-IVF™ PLUS in each step. Transfer the oocytes to the dishes with equilibrated G-IVF™ PLUS and return the dishes to the incubator immediately.

Insemination:

Transfer the oocytes to insemination dishes and leave at +37 °C and 6 % CO₂ over night.

Alternatively: Add equilibrated sperm suspension to equilibrated dishes with the oocytes already present. It is recommended to inseminate in 0.5–1 mL volume without oil overlay. If overlay is used, droplets of at least 100 µL are recommended.

Specifications

Sterile filtered	SAL 10 ³
Mouse Embryo Assay (1 cell)	
[% expanded blastocyst on day 5]	≥ 80
Bacterial endotoxins (LAL assay)	
[EU/mL]	< 0.25

LOT specific test results are available on the Certificate of Analysis provided with each delivery.

Precautions and warnings

Do not use G-IVF™ PLUS if it appears cloudy.

G-IVF™ PLUS contains human serum albumin and gentamicin.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested for antibodies to HIV, HBV, HCV and HTLV I/II

and non-reactive for HbsAg, HCV RNA and HIV-1 RNA and syphilis. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Vitrolife recommends that media be opened and used only with aseptic technique.

The risk of reproductive toxicity and developmental toxicity for IVF media, including Vitrolife's IVF media, have not been determined and are uncertain.

Not for injection.

Caution: Federal (US) law restricts this device to sale by or on the order of a physician.

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Vitrolife Sweden AB, Box 9080, SE-400 92 Göteborg, Sweden
Vitrolife Inc., South Inca Street, Englewood CO 80110, USA
Vitrolife Pty. Ltd., 26 Tyne Street, Carlton VIC 3053, Australia
Vitrolife Sweden AB